



Tc-99 radiochemistry in Risø/RAS

Jernström, Jussi

Publication date:
2009

[Link back to DTU Orbit](#)

Citation (APA):
Jernström, J. (Author). (2009). Tc-99 radiochemistry in Risø/RAS. Sound/Visual production (digital)

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.


Tc-99 RADIOCHEMISTRY IN RISØ/RAS

Jussi Jernström, Radiation Research Division, Risø DTU



Technetium as an element

- Main isotope Tc-99
 - Long half life; 211000 y
- Oxidation states from -1 to +7
- Almost totally anthropogenic element
- Occurs in nature only minute amounts by
 - spontaneous fission product in uranium ore
 - neutron capture in molybdenum ores
- Chemically between rhenium and manganese; chemistry close to Re



- Main emission of Tc-99 is β^- with 294 keV


Tc-99 in the environment

- Main sources of environmental Tc-99 are:
 - nuclear fuel cycle; primarily by fuel reprocessing (Sellafield and La Hague)
 - In Sellafield EARP technology started in 1994 increased Tc releases
 - La Hague releases decreased after the 1980's
 - nuclear weapon tests
 - diagnostic nuclear medicine (minor)
- Tc-99 is present mainly in marine environment due to controlled process discharges to seas
- Exists in oxic marine conditions as pertechnetate ion TcO_4^-
 - very soluble species; low K_d to sediments
- Reduces from Tc(VII) to Tc(IV) in reducing conditions
 - affinity to particles increases

Tc-99 in the marine environment

- Interest towards Tc-99 research is high
 - mobile Tc (TcO_4^-) can be used as marine tracer
 - risks of escape of mobile Tc from nuclear waste repositories
- Levels of Tc-99 in sea water
 - fallout background in Atlantic Ocean $\sim 0,005 \text{ Bq/m}^3$ [Dahlgaard et al. 1995]
 - North Sea and Norwegian coast (1996-1997): $0,8\text{-}8,5 \text{ Bq/m}^3$ [Brown et al. 1999]
 - Swedish west coast (2001) $\sim 1 \text{ Bq/m}^3$ [Lindahl et al. 2003]
 - Baltic Sea; Finnish coast (1999) $< 0,16 \text{ Bq/m}^3$ [Ikaheimonen et al. 2002]

Water circulation in the North Sea area



Anon 1993

Tc-99 speciation in the marine environment

- Behaviour of Tc in anoxic/reducing waters
 - Baltic Sea concentration drop
 - reduction of Tc and transfer from water column
- High Tc K_d values observed in some sediments (Techa river)
 - organic matter
 - anoxic conditions
- High uptake of Tc in brown algae
 - speciation
- New future activity in Risø/RAS

Tc-99 in soils


- Under oxidizing conditions in near-surface groundwater and soil Tc as TcO_4^-
 - soluble and stable over wide pH range
 - available for transfer to plants and biological processes
 - uptake mechanism?
- Insoluble species under reducing conditions
 - complexation with organic matter
 - some soil types together with wet conditions reduce Tc to be less available
 - soils with anion exchange capacity
 - anaerobic bacterial reduction
 - rice paddy fields

Monitoring of Tc-99 by Risø/RAS

- Long history of monitoring Tc
- Water and seaweed samples quarter-yearly from Danish waters
 - mainland
 - Faroe Islands
 - Greenland
- Water samples 200 L / analysis
- Seaweed samples
 - gram level for analysis



Tc-99 in seawater around Sjælland



Technetium-99 in seawater around Zealand, 1999-2009

Determination of Tc-99

- By radiometric methods using beta counting
 - liquid scintillation counting (LSC)
 - gas flow GM counting
 - neutron activation analysis
- By mass spectrometry
 - AMS
 - RIMS
 - ICP-MS

Radiometric methods vs. mass spectrometric methods

Table 9 – Comparison of radiometric and mass spectrometric methods for the determination of ^{99}Tc


Sample	Detection method	Target preparation	Detection limit	Sep. time ^a	Count time
Environmental sample	GM detector	Tc on disk	1.5 mBq	1–2 days/7 days	3–4 h
Environmental sample	LSC	Tc solution	17 mBq	1–2 days/7 days	2 h
Geological samples	RIMS	TcO ₂	20 μBq	1–2 days	30–60 min
Water	AMS	Tc in Al ₂ O ₃ or Nb ₂ O ₅	6–10 μBq	1 days	20–40 min
Sediment, seawater	ICP-SFMS	Tc solution	0.16–0.29 mBq	1–2 days	10–20 min
Seawater	ICP-MS with ETV	Tc in solution	0.18 mBq	1–2 days	20–40 min
Environmental sample	ICP-QMS	Tc in solution	10 mBq	1–2 days	20–40 min
Environmental sample	ICP-SFMS on-line	Tc in solution	0.05 mBq mL ⁻¹ or 0.2 mBq	1 h/4–5 h	10–20 min

Hou and Roos 2008

Tc-99 determination in water samples by beta counting

- High purification from other radionuclides needed
- Main analytical steps (Chen method)
 - purification and preconcentration by anion exchange
 - Tc absorbs in column in neutral or alkaline pH
 - washing away especially Ru-103 and Ru-106
 - purification by Fe-hydroxide precipitation
 - Tc does not precipitate
 - removal of U, Th, Pu, Am, Po, Np
 - further purification with solvent extraction
 - by using TIOA (= tri-iso-octylamine) in xylene
 - back-extraction with NaOH
 - electrodeposition of Tc on steel disk for GM counting


A procedure used for 200 L water samples in Risø/RAS



Tc-99 determination in seaweed samples by beta counting

- High purification from other radionuclides needed
- Main analytical steps (Chen method)
 - purification by Fe-hydroxide precipitation
 - Tc(VII) does not precipitate
 - removal of U, Th, Pu, Am, Po, Np
 - purification and preconcentration by anion exchange
 - Tc(VII) absorbs in column in neutral or alkaline pH
 - washing away especially Ru-103 and Ru-106
 - further purification with solvent extraction
 - by using TIOA in xylene
 - back-extraction with NaOH
 - electrodeposition of Tc on steel disk for GM counting


A procedure used for seaweed samples in Risø/NUK



Tc-99 determination in soil by mass spectrometry

- Main analytical steps:
 - Fe hydroxide precipitation
 - removal of impurities
 - Tc(VII) remains in solution
 - Anion exchange
 - purification and preconcentration
 - Tc(VII) absorbs in the resin
 - TEVA extraction chromatography
 - purification from Ru and Mo
 - Tc(VII) absorbs in resin
 - measurement with ICP-MS

A procedure for determining Tc-99 in water/soil with ICP-MS



- Procedure takes 1-2 working days
- Determination limit 0,3 mBq/m³
 - depends on ICP used

Interfering species in determination of Tc with ICP-MS

Table 8 – Potential interfering species on mass 99 in ICP-MS

Isobar/tailing	Oxide	Hydride	Argide	Chloride
⁹⁹ Ru (12.7%)	⁸³ Kr ¹⁶ O (11.5%)	⁹⁸ MoH (23.8%)	⁵⁹ Co ⁴⁰ Ar (100%)	⁶² Ni ³⁷ Cl (3.6%)
⁹⁸ Ru (1.1%)	⁸¹ Br ¹⁸ O (49.5%)	⁹⁸ RuH (1.9%)	⁶³ Cu ³⁶ Ar	⁶⁴ Zn ³⁵ Cl (48.6%)
¹⁰⁰ Ru (12.6%)	⁶⁷ Zn ¹⁶ O ₂ (4.1%)		⁴³ Ca ¹⁶ O ⁴⁰ Ar	
⁹⁸ Mo (24.1%)	⁵¹ V ¹⁶ O ₃ (99.8%)		⁴⁰ Ca ¹⁸ OH ⁴⁰ Ar	
¹⁰⁰ Mo (9.6%)				

- In low-Tc samples determination limit is hampered due to interfering species as the signal intensity is small
- Main interferences Ru-99 and hydride of Mo-98
 - Ru can be calculated from Ru-102
 - still the presence of Ru is a problem in low activity samples

Other methods for determination of Tc

- Lund method
 - wet ashing + ashing (with ammonia)
 - Fe hydroxide precipitation
 - removal of radiocesium with AMP
 - oxidation of Tc to Tc(VII) and extraction with TBP
 - HF to complex other elements
 - back-extraction with NaOH
 - electroplating of Tc
 - determination limit ~ 1,2 mBq/sample


Other methods for determination of Tc

- Southampton method
 - ashing (up to 550 °C) -> dissolution of Tc with 8 M HNO₃
 - Fe hydroxide precipitation
 - anion exchange
 - removal of calcium, Co-60, Ni-63, Sr-90, Cs-137
 - extraction with TnOA in xylene
 - Ru in water phase
 - organic phase with Tc to LSC solution
 - counting after 1 week (Tc-99m dies away)
 - determination limit 1,7 mBq/m³

Other methods for determination of Tc

- Harvey method
 - filtered sea water sample; addition of Re carrier
 - Re used as yield determinant
 - anion exchange
 - destruction of Tc-containing organic resin by heating
 - dissolution in HCl
 - Fe hydroxide precipitation
 - anion exchange
 - elution of Tc and Re with alkaline sodium perchlorate solution
 - sulphide precipitation of Tc and Re -> dissolution
 - precipitation of Tc and Re tetraphenylarsenium salt
 - gravimetric yield determination
 - measurement with proportional counter
 - determination limit for 100 L water sample 140 mBq/m³

Ion chromatography





- Why?
 - to elute Tc fraction being free from Mo and Ru
- Suitable complexing agent to create separation
 - suitable chemistry (pH and ionic strength)
 - suitable elution speed


Dionex CS5A column and its previous use for radionuclides



- Mixed-bed column for transition metals
- Separation of Am and Cm from lanthanides for alpha spectrometry
- Separation of Am and Pu




Purification of Tc from Mo and Ru




- Use of complexing agents
 - EDTA
 - Dipicolinic Acid (PDCA)
 - Oxalic Acid
- Chromatographic separation of Tc prior to MS
- Sample free from interfering Mo and Ru
- Tc both in reduced and oxidized forms

Complexing agents


EDTA



Oxalic acid




PDCA




Reduction of Tc

- 20 mM ascorbic acid + 20 mM hydroxylamine
- 25 min reduction time
- EDTA and oxalic acid positive
- PDCA does not elute all Tc


[EDTA] 0,1 M; pH 7,0; [Na₂HPO₄] 50 mM; [Na₃PO₄] 50 mM
[Ascorbic acid] 20 mM; [Hydroxylamine] 20 mM



[Oxalic acid] 0,15 M; pH 4,2; [CH₃COOH] 75 mM; [CH₃COONa] 75 mM
[Ascorbic acid] 20 mM; [Hydroxylamine] 20 mM







[PDCA] 50 mM; pH 4,65; [CH₃COOH] 75 mM; [CH₃COONa] 75 mM
[Ascorbic acid] 20 mM; [Hydroxylamine] 20 mM




No reducing agents; Tc as Tc(VII)


- EDTA and oxalate show no ability for elution of Tc or Mo and Ru
- PDCA complexes strongly -> full separation from Ru
 - No elution of Mo



Results with soil sample



- elution with oxalate



- elution with PDCA

Thank you for your attention!