



## Use of ICP-MS in analysing radioisotopes

**Roos, Per**

*Publication date:*  
2009

[Link back to DTU Orbit](#)

*Citation (APA):*  
Roos, P. (Author). (2009). Use of ICP-MS in analysing radioisotopes. Sound/Visual production (digital)

---

### General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

# Use of ICP-MS in analysing radioisotopes

Per Roos

Risø National Laboratory for Sustainable Energy, Technical University of Denmark





Fig. 1. Overview of application fields for determination of long-lived radionuclides.

# Inductively Coupled Plasma – Mass Spectrometry (ICP-MS)

## History




- ICP-AES in 1960'ies
- Plasma & MS merging in 1974 using a DC-plasma, Surrey UK.  
(A.L. Gray, A plasma source for mass analysis, *Proc. Soc. Anal. Chem*, 11, 182-183).
- ICP & MS merging in 1980, Ames Lab, USA  
(R.S. Houk et.al., Inductively coupled argon plasma as an ion source for mass spectrometric determination of trace elements. *Anal. Chem*, 52, 2283-2289.).
- Microwave plasma – MS, Toronto Canada.  
(D.J. Douglas & J.B. French, Elemental analysis with a microwave induced plasma quadrupole mass spectrometer system. *Anal. Chem.* 53, 37-41.)
- First commercial instruments in 1983.  
PlasmaQuad by VG Isotopes, "UK-system"  
Elan by Sciex, "Canada system"


## Inductively Coupled Plasma

- Electrodeless discharge formed in a flowing gas stream enclosed in a quartz torch. Sustained by the induction of electrical energy using a coil around the torch. Atmospheric pressure!
- "Plasma": Langemuir in 1928 adopted the word to describe a discharge through a gas confined at low pressure.
- "Any volume of an ionized gas where the total charge is zero"

## Lightning the plasma

- Seed electrons to the flowing gas (piezoelectric spark)
- Apply a high-frequency electromagnetic field (20-50 MHz)
- Electrons will be accelerated back- and forwards in oscillating RF-field
- Multiplication of free electrons – plasma reached in a few milliseconds.





**Figure 1.4** (A) Schematic diagram of a quadrupole ICPMS, the Elan 6000. (Courtesy of the Perkin-Elmer/Sciex Corporation.) (B) Schematic diagram of a magnetic sector ICPMS, the JMS-PLASMAX2. (Courtesy of JOEL Incorporated.)


## Risø ICP-MS instruments

Thermo X-series II  
Quadrupole ICP-MS



Plasma Trace 2  
Sector field ICP-MS

## Torch, loadcoil and interface







# Torch & load coil

- Three concentric coaxial glass/quartz tubes
- Outer tube: About 25mm longer than the two others.
- Cool gas: Rotating 12-16 L min<sup>-1</sup>
- "Auxiliary" gas tube (plasma gas): 0.6-1.4 L min<sup>-1</sup>
- Injector or nebuliser gas: 0.5-1 L min<sup>-1</sup>
- Induction coil (load coil);
- Water cooled
- Copper
- Part of the oscillatory circuit called the RF-generator.
- 20-40 MHz
- Power output 1-2.5 kW.












## The plasma – ms interface

- Sampling and skimming process
- In early days only sampler- 0.05mm
- Sample orifice must be large enough to sample centre of plasma while at the same time small enough not to disturb the plasma.






**Table 1.3** Distribution of ionization energies among the elements for singly and doubly charged ions. Grouped in steps of 1 eV.

Ionization Energy (eV)	Elements
< 7	Li, Na, Al, K, Ca, Sc, Ti, V, Cr, Ga, Rb, Sr, Y, Zr, Nb, In, Cs, Ba, La, Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu, Hf, Tl, Ra, Ac, Th, U
7–8	Mg, Mn, Fe, Co, Ni, Cu, Ge, Mo, Tc, Ru, Rh, Ag, Sn, Sb, Ta, W, Re, Pb, Bi
8–9	B, Si, Pd, Cd, Os, Ir, Pt, Po
9–10	Be, Zn, As, Se, Te, Au
10–11	P, S, I, Hg, Rn
11–12	C, Br
12–13	Xe
13–14	H, O, Cl, Kr
14–15	N
15–16	Ar
> 16	He, F, Ne

*2<sup>+</sup> Ions*

Ba, Ce, Pr, Nd, Ra  
Ca, Sr, La, Sm, Eu, Tb, Dy, Ho, Er  
Sc, Y, Gd, Tm, Yb, Th, U, Ac  
Ti, Zr, Lu  
V, Nb, Hf  
Mg, Mn, Ge, Pb  
All other elements



**Figure 1.7** Degree of ionization  $\alpha$  versus ionization energy for singly charged ions in the ICP. Calculated for a representative selection of 59 elements from the Saha equation, assuming values of  $T_i$  of 8000 K and  $n_e$  of  $2.5 \times 10^{15} \text{ cm}^{-3}$ .

H 0.1																	He
Li 100	Be 75											B 58	C 5	N 0.1	O 0.1	F 9e-4	Ne 6e-8
Na 100	Mg 98											Al 98	Si 85	P 33	S 14	Cl 0.9	Ar 0.04
K 100	Ca 99	Sc 100	Ti 99	V 99	Cr 98	Mn 95	Fe 96	Co 93	Ni 91	Cu 90	Zn 75	Ga 98	Ge 90	As 52	Se 33	Br 5	Kr 0.6
Rb 100	Sr 96	Y 98	Zr 99	Nb 98	Mo 98	Tc	Ru 96	Rh 94	Pd 93	Ag 93	Cd 85	In 99	Sn 96	Sb 78	Te 66	I 29	Xe 8.5
Cs 100	Ba 91	La 90	Hf 98	Ta 95	W 94	Re 93	Os 78	Ir	Pt 62	Au 51	Hg 38	Tl 100	Pb 97	Bi 92	Po	At	Rn
Fr	Ra	Ac															
			Ce 98	Pr 90	Nd 99	Pm	Sm 97	Eu 100	Gd 93	Tb 99	Dy 100	Ho	Er 99	Tm 91	Yb 92	Lu	
			Th 100	Pa	U 100	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lw	

**Figure 8.1** Calculated degrees of ionization,  $[M^+]/([M^+] + [M])$ , expressed as percentages. An ionization temperature of 7500 K and electron number density of  $1 \times 10^{15} \text{ cm}^{-3}$  were used. (From Reference 10, with permission.)

Table 7

Detection limits for actinides in high-purity water by double focusing sector fields ICP-MS

Nuclide	Half life [years]	“ELEMENT” [129]		“Plasma Trace” [175] [pg L <sup>-1</sup> ]
		[pg L <sup>-1</sup> ]	[kBq L <sup>-1</sup> ]	
<sup>230</sup> Th	$7.5 \times 10^4$	0.08	$6.2 \times 10^{-8}$	20
<sup>232</sup> Th	$1.4 \times 10^{10}$	0.1	$4.1 \times 10^{-13}$	
<sup>233</sup> U	$1.6 \times 10^5$	0.07	$2.5 \times 10^{-8}$	
<sup>238</sup> U	$4.5 \times 10^9$	0.2	$2.4 \times 10^{-12}$	20
<sup>237</sup> Np	$2.1 \times 10^{16}$	0.05	$1.3 \times 10^{-9}$	2
<sup>239</sup> Pu	$2.4 \times 10^4$	0.04	$9.4 \times 10^{-8}$	5
<sup>243</sup> Am	$7.9 \times 10^3$	0.05	$3.6 \times 10^{-7}$	

**Table 9.15** The detection limits (DL) obtained using magnetic sector-ICP-MS with ultrasonic nebuliser

Nuclide	Half-life (years)	CPS per pg g <sup>-1</sup>	DL (fg g <sup>-1</sup> )	DL (Bq g <sup>-1</sup> )
<sup>99</sup> Tc	$2.14 \times 10^5$	$2.28 \times 10^2$	2.6	$1.60 \times 10^{-6}$
<sup>226</sup> Ra	$1.60 \times 10^3$	$1.01 \times 10^4$	0.1	$3.70 \times 10^{-6}$
<sup>232</sup> Th	$1.41 \times 10^{10}$	$7.78 \times 10^3$	0.2	$8.00 \times 10^{-13}$
<sup>237</sup> Np	$2.14 \times 10^6$	$5.83 \times 10^3$	0.2	$5.20 \times 10^{-9}$
<sup>238</sup> U	$4.47 \times 10^{10}$	$8.02 \times 10^3$	0.2	$2.40 \times 10^{-12}$
<sup>239</sup> Pu	$2.41 \times 10^4$	$6.21 \times 10^3$	0.2	$4.28 \times 10^{-7}$
<sup>240</sup> Pu	$6.56 \times 10^3$	$6.21 \times 10^3$	0.2	$1.57 \times 10^{-6}$
<sup>243</sup> Am	$7.37 \times 10^3$	$6.56 \times 10^3$	0.2	$1.40 \times 10^{-6}$



**Table 5.** Application of ICP-MS in trace and isotope analysis of long-lived radionuclides in nuclear and environmental materials, radioactive waste and medical samples

Samples	Method	Radionuclides	Detection limits/isotope ratio	References
Fission product solution	ETV-ICP-MS	$^{79}\text{Se}$	13 ng/l	Compte et al. [72]
Thermal water	ICP-QMS + USN (PlasmaQuad PQ2) extraction, ion exchange	$^{226}\text{Ra}$	2 fg/g	Joannon and Pin [69]
Mineral water	ICP-SFMS after separation	$^{226}\text{Ra}$	2 fg/g	Zoriy et al. [73]
Radioactive waste solution	ICP-SFMS ELEMENT, (USN)	$^{226}\text{Ra}$ , $^{230}\text{Th}$ , $^{233}\text{U}$ , $^{237}\text{Np}$ , $^{239}\text{Pu}$ , $^{241}\text{Am}$	0.05 pg/l ( $^{241}\text{Am}$ ) 0.04 pg/l ( $^{239}\text{Pu}$ ) $^{236}\text{U}/^{238}\text{U} \sim 0.02$ (flow injection)	Becker and Dietze [7]
Radioactive waste solution	ICP-SFMS ELEMENT, (DIHEN)	$^{226}\text{Ra}$ , $^{230}\text{Th}$ , $^{230}\text{Th}$ , $^{233}\text{U}$ , $^{237}\text{Np}$ , $^{238}\text{U}$ , $^{241}\text{Am}$	0.1 pg/l ( $^{241}\text{Am}$ ) 0.1 pg/l ( $^{237}\text{Np}$ ) $^{236}\text{U}/^{238}\text{U} \sim 0.00002-0.0001$	McLean et al. [30]
Urine	ICP-SFMS (ELEMENT) after separation (1 l urine)	$^{239}\text{Pu}$ , $^{240}\text{Pu}$	1 fg/l $^{240}\text{Pu}/^{239}\text{Pu} \sim 0.14$ (synth. sample)	Zoriy et al. [54]
Urine	ICP-SFMS (ELEMENT)	$^{90}\text{Sr}$	0.4 pg/l (ELEMENT)	Vonderheide et al. [53]
	ICP-CC-MS (Platform)		2 ng/l (Platform)	
Ground water	ICP-SFMS (ELEMENT)	$^{90}\text{Sr}$ , $^{239}\text{Pu}$ , $^{240}\text{Pu}$	11 pg/l ( $^{90}\text{Sr}$ ); 0.1 pg/l (Pu); $^{240}\text{Pu}/^{239}\text{Pu} \sim 0.17$	Zoriy [71]
Sea water	ICP-SFMS (ELEMENT)	$^{239}\text{Pu}$ , $^{240}\text{Pu}$	$\approx 0.1$ fg/l (ELEMENT) $\approx 0.01$ fg/l (NuPlasma) $^{240}\text{Pu}/^{239}\text{Pu} \sim 0.17$	Becker et al. [52]

## Radioisotopes – in general focus on:

- Sensitivity (transmission of ions to detector)
- Polyatomic interferences
- Abundance sensitivity
- Stability (isotope ratios)

# Improving sensitivity

- Matrix removal from sample (separation chemistry)
- High efficiency sample introduction systems
- Instrument optimisation

## Pu-analysis at sub-mBq to $\mu$ Bq levels


- Instrument background.
- Procedure blank (Pu)
- Contribution from tracer.
- Memory effects.
- Polyatomic species (including the  $\text{UH}^+$ )
- Abundance sensitivity.

## Transport efficiency

- 1 mBq ml<sup>-1</sup> of <sup>242</sup>Pu ( $1.7 \cdot 10^{10}$  atoms ml<sup>-1</sup>)
  - With USN and 0.4 ml min<sup>-1</sup> uptake rate produces a count rate of about 10 000 cps.
- 
- $\Rightarrow$  **1 atom per 12 500 detected !**

## Efficiency alpha vs. ICP-MS

- **Alpha**: Geometric efficiency about 25%
- **ICP-MS**: Efficiency about 1: 10 000
- **Alpha counting** of  $^{239}\text{Pu}$  (24 400 y). Two months counting  $\Rightarrow$  about  $10^{-6}$  of sample atoms have decayed (for  $^{238}\text{U}$  about  $10^{-12}$ ).
- **ICP-MS** measurement of  $^{239}\text{Pu}$  takes about 10 minutes.  $10^{-4}$  of sample atoms are detected.



# Ways to improve sensitivity

## Shielded torch

Fig. 15. Plasma shielded torch (GE—guard electrode, courtesy of Finnigan MAT).






Fig. 16. Comparison of sensitivity of long-lived radionuclides in DF-ICP-MS with and without shielded torch. (diagonal lines, with shielded torch, vertical lines, without shielded torch).

# High efficiency nebulisers

(B)








## Sensitivity & sample introduction system

	Concentric nebuliser + chilled cyclonic spraychamber	USN
cps per pg s <sup>-1</sup> to plasma (uranium)	1.2*10 <sup>6</sup>	1.1 *10 <sup>8</sup>



# Electrothermal vaporization (ETV)

.....also for speciation studies






**Figure 9.7** Diagram of graphite tube ETV-ICPMS atomizer.




**Figure 9.8** Signal responses for typical analytes by ETV-ICPMS. (From Reference 91, with permission.)

## ....APEX-system



The importance of removing uranium in Pu-measurements (hydride generation and spillover)



# The UH<sup>+</sup> and ways to minimize it




Fig. 2. Uranium hydride ion production in ICP-MS using PFA 100 nebulizer (a) without and (b) with the APEX system.

UH<sup>+</sup>/U

Wet plasma (Cyclonic spraychamber)	Dry plasma (USN)
63 ppm	23 ppm

Hydrides important to minimize for  $^{236}\text{U}$  analysis.






Fig. 4.  $^{236}\text{U}/^{238}\text{U}$  vs.  $^{235}\text{U}/^{238}\text{U}$  ratios in soil samples collected in different regions of Belarus.



# Direct Injection High Efficiency Nebuliser (DIHEN)




**Figure 3.22** Schematic diagram of the DIHEN interfaced with a demountable plasma torch. (Revised from Reference 269, with permission.)


# Best use of a limited amount of sample

## Efficiency versus uptake rate

(Concentric nebuliser - chilled cyclonic spraychamber)




**Fraction of sample reaching torch**  
(Concentric nebuliser - chilled cyclonic spray chamber)



# Efficiency of sample reaching plasma

(Concentric nebuliser - chilled cyclonic spray chamber)



# Isotope ratio measurements

- Mass discrimination
- Dead time corrections
- Reference materials !

Table 2


Precision and accuracy of DIHEN-ICP-SFMS measuring uranium isotope standard solutions NIST U005, U350, U930 and U(nat) ( $U_{\text{conc}} = 0.1$  ppb,  $n = 6$ )

	Isotope ratios		Precision (R.S.D., %)	Accuracy (%)
	Measured	Certified		
U005				
$^{235}\text{U}/^{238}\text{U}$	0.0049435	0.0049194	1.4	0.5
$^{234}\text{U}/^{238}\text{U}$	0.0000218	0.0000219	1.8	−0.12
$^{236}\text{U}/^{238}\text{U}$	0.0000472	0.0000468	4.3	0.9
U350				
$^{235}\text{U}/^{238}\text{U}$	0.540114	0.546488	1.6	−1.2
$^{234}\text{U}/^{238}\text{U}$	0.003863	0.003879	3.2	−0.41
$^{236}\text{U}/^{238}\text{U}$	0.002644	0.002598	1.3	1.8
U930				
$^{235}\text{U}/^{238}\text{U}$	17.3472	17.3487	2.0	−0.01
$^{234}\text{U}/^{238}\text{U}$	0.2005	0.2010	1.9	−0.23
$^{236}\text{U}/^{238}\text{U}$	0.03783	0.03768	2.8	0.4
U(nat)				
$^{235}\text{U}/^{238}\text{U}$	0.007227	0.007253	1.1	−0.3
$^{234}\text{U}/^{238}\text{U}$	0.000054	0.000055	0.1	−0.52

## Problems in isotope ratio measurements (single collectors)


- Sample introduction system
- Plasma instability
- Dead time corrections
- Mass fractionation
- Interferences
- Abundance sensitivity (peak tailing)

## Mass discrimination at interface



**Plate 2** Under high charge density conditions, the Debye spheres of the charged particles overlap and the particles experience modified fields (they are “shielded” from the applied electrostatic field). Because of their higher mobility, electrons preferentially diffuse to the walls where an electron sheath forms. Due to the electron diffusion, a net positive charge density remains on axis and the resultant space charge field is positive. See color plate.

# Mass discrimination



**Fig.2** Mass discrimination per mass unit in ICP-QMS. —■— Platform ICP-MS; —▲— Elan 5000 [1] (rf power 1350 W, nebulizer gas 0.86 L min<sup>-1</sup>; collision gas: He – 4 ml min<sup>-1</sup>; H<sub>2</sub> – 1.5 min min<sup>-1</sup>, hexapole bias 0 V)<sup>1</sup>






Fig. 5. Plutonium isotope ratio analysis at femtogram to nanogram levels by MC-ICP-CC-MS vs. TIMS (from Ref. [9] with permission).

**Table 2** Isotope ratio measurements by different ICP-MS methods

Isotope ratio	Method	Precision	Sample size	Reference
$^{235}\text{U}/^{238}\text{U}\sim 1$	ICP-QMS	0.028% <sup>a</sup>	255 pg	Platzner et al. [24]
$^{207}\text{Pb}/^{206}\text{Pb}\sim 0.9$	DRC-ICP-QMS	0.057% 0.027% <sup>a</sup>	440 ng 5.1 µg	Bandura and Tanner [44]
$^{235}\text{U}/^{238}\text{U}\sim 1$	DF-ICP-MS	0.026% 0.2%	255 pg 0.5 pg	Becker and Dietze [28]
$^{34}\text{S}/^{32}\text{S}\sim 0.05$	Hex-ICP-QMS	0.03%	50–250 µg	Mason et al. [40]
$^{137}\text{Ba}/^{138}\text{Ba}$	TOF-ICP-QMS	0.05%	250 ng	Vanhaecke et al. [34]
$^{208}\text{Pb}/^{206}\text{Pb}$	MC-ICP-MS	0.006%	100 mg	Rehkämper and Halliday [52]
$^{182}\text{W}/^{183}\text{W}\sim 2$	MC-ICP-MS	0.006%	0.5–1 µg	Lee and Halliday [50]
$^{63}\text{Cu}/^{65}\text{Cu}$	MC-ICP-MS	0.002%	100 ng	Maréchal et al. [54]

<sup>a</sup>long term precision

# Different separation chemistry relative to "alpha spectrometry"


Table 1  
Possible interferences for Pu isotopes and required mass resolution on ICP-SFMS

Nuclide	Molecular ions	Required mass resolution, $m/\Delta m$
$^{238}\text{Pu}$	$^{238}\text{U}^+$	193665
	$^{206}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}_2^+$	3874
	$^{208}\text{Pb}^{16}\text{O}^1\text{H}_2^{12}\text{C}^+$	3818
	$^{208}\text{Pb}^{14}\text{N}_2^1\text{H}_2^+$	4654
$^{239}\text{Pu}$	$^{238}\text{U}^1\text{H}^+$	36885
	$^{207}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}_2^+$	3817
	$^{208}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}^+$	3430
$^{240}\text{Pu}$	$^{238}\text{U}^1\text{H}_2^+$	19116
	$^{208}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}_2^+$	3774
$^{241}\text{Pu}$	$^{207}\text{Pb}^{16}\text{O}_2^1\text{H}_2^+$	3193
$^{242}\text{Pu}$	$^{208}\text{Pb}^{16}\text{O}_2^1\text{H}_2^+$	3159
$^{244}\text{Pu}$	$^{206}\text{Pb}^{12}\text{C}_3^1\text{H}_2^+$	3293
	$^{207}\text{Pb}^{12}\text{O}_3^1\text{H}^+$	3031


Table 1  
Possible interferences for  $^{90}\text{Sr}$  and Pu isotopes and required mass resolution on ICP-SFMS

Nuclide	Molecular ions	Required mass resolution ( $m/\Delta m$ )
$^{90}\text{Sr}$	$^{180}\text{W}^{2+}$	1370
	$^{180}\text{Hf}^{2+}$	1372
	$^{58}\text{Ni}^{16}\text{O}_2^+$	2315
	$^{74}\text{Ge}^{16}\text{O}^+$	10765
	$^{52}\text{Cr}^{38}\text{Ar}^+$	19987
	$^{50}\text{V}^{40}\text{Ar}^+$	49894
	$^{54}\text{Fe}^{36}\text{Ar}^+$	155548
	$^{50}\text{Ti}^{40}\text{Ar}^+$	158287
$^{239}\text{Pu}$	$^{90}\text{Zr}^+$	29877
	$^{238}\text{U}^1\text{H}^+$	36885
	$^{207}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}_2^+$	3817
	$^{208}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}^+$	3430
$^{240}\text{Pu}$	$^{238}\text{U}^1\text{H}_2^+$	19116
	$^{208}\text{Pb}^{16}\text{O}^{14}\text{N}^1\text{H}_2^+$	3774

# 2 ppm Pb + 5% HCl




# Tuning for Pu-analysis – we use uranium




# The need to remove uranium in Pu-analysis


## Peak tailing & peak shape



$^{235}\text{U}/^{238}\text{U}$  atom ratio vs.  $^{235}\text{U}$  counts




# Abundance sensitivity thorium, natural soil sample without radiochemistry







## Abundance sensitivity natural uranium in soil. No radiochemistry



Pu3-IAEA-384-7\_1




# Polyatomic interferences – Pu analysis (stripped steel disc)




# Uranium isotopes


- $^{238}\text{U}$ ,  $^{235}\text{U}$ ,  $^{234}\text{U}$  as natural U-isotopes.
- $^{236}\text{U}$  ( $^{233}\text{U}$ ) as artificial U-isotopes.
- Complicated isotope signature in recycled uranium require MS.




## Depleted and recycled uranium



# $^{238}\text{U}$ raw counts




## Depleted and recycled uranium – alpha spectrometry






## Pu without tracer – alpha spectrometry



# Pu-isotopes Irish Sea Sediments - ICP-MS



Pu-990424\_1

ICPMS-10  
Pu-030999

## 50 L seawater, Baltic (Møen), Surface




Table 7  
 Characteristics of different ICP-MS instruments using microcentric nebulizer with desolvator (Aridus, Cetac Technologies)

	Absolute sensitivity (counts/atom)	$\text{UH}^+/\text{U}^+$	Abundance sensitivity
ICP-QMS <sup>a</sup>	$4 \times 10^{-5}$	$3 \times 10^{-6}$	$6 \times 10^{-7}$
ICP-SFMS	$7 \times 10^{-4}$	$3 \times 10^{-5}$	$5 \times 10^{-6}$
MC-ICP-MS <sup>b</sup>	$5 \times 10^{-4}$	$5 \times 10^{-5}$	$3 \times 10^{-7}$

<sup>a</sup> ELAN 6000.

<sup>b</sup> Halicz [94].

## Determinations of $^{240}\text{Pu}/^{239}\text{Pu}$ atom ratio in Normass reference soil

(Leached soil solution)

	$^{240}\text{Pu}/^{239}\text{Pu}$
Lab 1	$0.187 \pm 0.004$
Lab 2	$0.193 \pm 0.004$
Lab 3	$0.187 \pm 0.006$
Kelley et.al. 'Roskilde soil' (Sci.Tot. Env. 237/238, 1999, 483-500)	$0.1904 \pm 0.0012$


<b>Run No.</b>	<b>Atom <math>^{240}\text{Pu}/^{239}\text{Pu}</math></b>	<b>Atom <math>^{237}\text{Np}/^{239}\text{Pu}</math></b>	<b><math>^{239+240}\text{Pu}</math> mBq/g ICP</b>	<b><math>^{239+240}\text{Pu}</math> mBq/g alfa</b>
<b>Soil-1</b>	<b>0.1908</b>	<b>0.283</b>	<b>0.240</b>	<b>0.251</b>
<b>Soil-2</b>	<b>0.1818</b>	<b>0.362</b>	<b>0.246</b>	<b>0.239</b>
<b>Soil-3</b>	<b>0.1969</b>	<b>0.329</b>	<b>0.239</b>	<b>0.238</b>
<b>Soil-4</b>	<b>0.1965</b>	<b>0.318</b>	<b>0.238</b>	<b>0.235</b>
<b>Soil-5</b>	<b>0.1895</b>	<b>0.294</b>	<b>0.227</b>	<b>0.231</b>
<b>Soil-6</b>	<b>0.1887</b>	<b>0.304</b>	<b>0.239</b>	<b>0.238</b>
<b>Soil-7</b>	<b>0.1928</b>	<b>0.316</b>	<b>0.240</b>	<b>0.246</b>
<b>Soil-8</b>	<b>0.1790</b>	<b>0.371</b>	<b>0.240</b>	<b>0.233</b>
			<b>0.238±0.0018</b>	<b>0.238±0.0087</b>
<b>Atom ratio:</b>	<b>0.1895±0.0023</b>	<b>0.322±0.01</b>		
<b>Roskilde (Kelley et.al.)</b>	<b>0.1904±0.0012</b>	<b>0.531±0.013</b>		

Run No.	$^{238}\text{U}$	$^{239}\text{Pu}$	$^{240}\text{Pu}$	Atom $^{240}\text{Pu}/^{239}\text{Pu}$
1	1223	5000	911	0.1822
2	1180	5310	990	0.1864
3	1143	5368	923	0.1719
4	1082	5196	955	0.1838
5	1135	5111	964	0.1886
6	1058	5185	957	0.1846
7	1069	5140	979	0.1905
8	1117	5174	938	0.1813
9	980	4937	975	0.1975
10	1081	5004	967	0.1932
11	1031	5013	959	0.1913
12	991	5080	968	0.1906
13	1012	4985	894	0.1793
14	958	5179	966	0.1865
15	1020	4899	948	0.1935
Atom ratio:				0.1868±0.0065 (1 std.dev)

# Collision- and reaction cells

## CCT<sup>ED</sup> - Kinetic Energy Discrimination

Collisional retardation / energy filtering





## Collision cell and I-129 analysis






Figure 9.15 Cell gas optimisation for  $^{129}\text{Xe}$  background.



**Figure 1.** Mass spectra for samples containing 1 ppb single-element standards of U, Pu, and Am (see text for impurities in Am sample), in standard mode without reaction gas.



**Figure 3.** Mass spectra for samples containing 1 ppb single-element standards of U, Pu, and Am with carbon dioxide gas (0.5 Ar-equivalent sccm, 0.25 true sccm).