



## A peristaltic pump driven $^{89}\text{Zr}$ separation module

**Siikanen, J.; Peterson, M.; Tran, T.; Roos, Per; Ohlsson, T.; Sandell, A.**

*Published in:*  
A I P Conference Proceedings Series

*Publication date:*  
2012

*Document Version*  
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

*Citation (APA):*  
Siikanen, J., Peterson, M., Tran, T., Roos, P., Ohlsson, T., & Sandell, A. (2012). A peristaltic pump driven  $^{89}\text{Zr}$  separation module. *A I P Conference Proceedings Series*, 1509, 206-210.

---

### General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

## **A peristaltic pump driven 89 Zr separation module**

J. Siikanen, M. Peterson, T. A. Tran, P. Roos, T. Ohlsson, and A. Sandell

Citation: [AIP Conference Proceedings](#) **1509**, 206 (2012); doi: 10.1063/1.4773969

View online: <http://dx.doi.org/10.1063/1.4773969>

View Table of Contents: <http://scitation.aip.org/content/aip/proceeding/aipcp/1509?ver=pdfcov>

Published by the [AIP Publishing](#)

---

### **Articles you may be interested in**

[A valve-less microfluidic peristaltic pumping method](#)

Biomicrofluidics **9**, 014118 (2015); 10.1063/1.4907982

[The influence of suspended drops on peristaltic pumping](#)

Phys. Fluids **26**, 061902 (2014); 10.1063/1.4882263

[Peristaltic pumping and irreversibility of a Stokesian viscoelastic fluid](#)

Phys. Fluids **20**, 073101 (2008); 10.1063/1.2963530

[Peristaltically driven channel flows with applications toward micromixing](#)

Phys. Fluids **13**, 1837 (2001); 10.1063/1.1377616

[Miniature peristaltic vacuum pump for use in portable instruments](#)

J. Vac. Sci. Technol. A **14**, 2610 (1996); 10.1116/1.579988

---

# A Peristaltic Pump Driven $^{89}\text{Zr}$ Separation Module

J. Siikanen<sup>a,b</sup>, M. Peterson<sup>b</sup>, T.A. Tran<sup>a,c</sup>, P. Roos<sup>d</sup>, T. Ohlsson<sup>b</sup> and A. Sandell<sup>b</sup>

<sup>a</sup>*Department of Medical Radiation Physics, Clinical Sciences Lund, Skånes University Hospital Lund, Barngatan 2:1, 221 85 Lund, Sweden*

<sup>b</sup>*Skånes University Hospital, Radiation Physics, Klinikgatan 7, 221 85 Lund, Sweden*

<sup>c</sup>*Lund University Bioimaging Center, Lund University, Faculty of Medicine, BMC D11221 84 Lund, Sweden*

<sup>d</sup>*DTU Nutech, Center for Nuclear Technologies, Nutech Radioecology and Tracer Studies, Technical University of Denmark, Frederiksborgvej 399 Dk-4000, Roskilde, Denmark*

**Abstract.** To facilitate the separation of  $^{89}\text{Zr}$  produced in yttrium foils, an automated separation module was designed and assembled. The module separates more than 85% of produced  $^{89}\text{Zr}$  – activity in 3 g foils in less than 90 min. About 10 % remains in the dissolving vial. The quality of the separated  $^{89}\text{Zr}$  activity was investigated for labeling of the HER2-binding monoclonal antibody fragment, trastuzumab-Fab.

**Keywords:**  $^{89}\text{Zr}$ , PET, Separation Module, Automation.

**PACS:** 28.60.+s, 87.57.uk, 87.57.un

## INTRODUCTION

$^{89}\text{Zr}$  ( $\beta^+=22.7\%$ ,  $t_{1/2}=3.27\text{ d}$ ) is a useful radiometal for immuno-PET, the *in vivo* imaging and quantification of monoclonal antibodies (mAbs) with positron emission tomography (PET). This is because  $^{89}\text{Zr}$  has relatively high  $\beta^+$  intensity and a physical half-life suitably matched to the pharmacokinetics of mAbs *in vivo* (days). Since the first utilization of  $^{89}\text{Zr}$  for labeling of intact antibodies was presented,<sup>1</sup>  $^{89}\text{Zr}$  has been given a great deal of attention.<sup>2</sup>

$^{89}\text{Zr}$  is accessible with hospital cyclotrons through the  $^{89}\text{Y}(\text{p,n})^{89}\text{Zr}$  reaction on natural yttrium ( $^{89}\text{Y}=100\%$ ). One common procedure to separate  $^{89}\text{Zr}$  from irradiated yttrium material is to use a solid-phase hydroxamate resin<sup>3,4</sup> from which  $^{89}\text{Zr}$  is eluted with oxalic acid as  $[^{89}\text{Zr}]\text{Zr-oxalate}$ . In 2010 Vosjan et al.<sup>5</sup> published a protocol describing the step by step procedure, where separated  $[^{89}\text{Zr}]\text{Zr-oxalate}$  is directly used to label mAbs or other proteins pre-conjugated with the commercially available bifunctional chelate p-isothiocyanatobenzyl-desferrioxamine (Df-Bz-NCS).

This work focuses on the separation of  $^{89}\text{Zr}$  activity from natural yttrium target material with a hydroxamate resin. For this purpose, an automated module was designed and assembled. An automated separation module is one important parameter in the streamlining of regular  $^{89}\text{Zr}$ -production. Additionally, personnel doses are decreased with an automated module.


## MATERIAL AND METHODS

### Equipment & Materials


The module consists of a peristaltic pump (Welco, 6-12 VDC) and 6 dual pinch valves (Takasago, 12 VDC) syringes and hydroxamate resin, which are mounted on an aluminum plate (Figure 1). On the upper side of the pinch valves, the syringes and the hydroxamate resin are connected with Luer Lock tips to a Pharmed BPT tube (1 mm i.d and 3 mm o.d). On the lower side of the pinch valves the tubes are connected to a common line with 3 way stopcocks (only used as tee-unions).

Pinch valves and pump are connected via an on/off 8 channel USB relay card (Quancom-USBREL8/LC) to 6-12 VDC lines. A PC lap-top with a LabVIEW (National Instruments) program is used to control the relay card. The separation program starts a sequence which consists of several steps. Each step has a user specified valve/pump setting and duration (Figure 2).

To decrease metal contamination low metal level HCl (Merck 30 % Ultrapur®) and water (Sigma Aldrich, Traceselect®) is used. MeCN is pretreated with chelex 100 resin (Bio-Rad Laboratories) to remove metal ions.



**FIGURE 1.** (a). The module loaded with resin, liquids and an irradiated foil in the dissolving vial.  
(b). Schematic overview of the set up.



**FIGURE 2.** LabVIEW Valve/Pump control-program

## Irradiation

$^{89}\text{Zr}$  is normally produced in  $\sim 50 \times 20 \times 0.64$  mm,  $\sim 3$ g, yttrium foils (Alfa Aesar 99.9%) with 45  $\mu\text{A}$ ,  $\sim 12.8$  MeV protons with a home-made solid multi-purpose target system.<sup>6</sup> One hour irradiation results in  $2150 \pm 197$  MBq ( $n=13$ ).

## Separation

The automated separation in this work is performed based on the separation described by Holland et al.<sup>7</sup> Briefly, a home manufactured 100 mg hydroxamate resin is manually activated, with MeCN, water, 2 M HCl and connected to the module. Also water and HCl-syringes are connected to the module before an irradiated yttrium foil is placed in the dissolving vial. The peristaltic pump is placed so that liquids on the right side of the pump can go to vials on the left side of the pump and vice versa by reversing the direction of the pump rotation i.e. switch the polarity of the motor. The program is started where 6 M HCl is added to dissolve the foil. Because of approximately ten times larger target masses in this work compared to Holland et al.,<sup>7</sup> water and HCl volumes are increased proportionally. Due to the violent reaction with HCl, foils are dissolved in several steps by pumping 1 ml HCl to the dissolving vial and waiting 75 sec before next addition (total of  $\sim 18$  ml). To collect fumes the exhaust is operating during the entire dissolving step by circulating water through a water jet vacuum pump. In the next step, water is added to ensure a  $< 2$  M HCl concentration. The solution is transferred to the resin where  $^{89}\text{Zr}$  is trapped. The resin is then washed with 2 M HCl and water to remove Y(III) ions and other impurities. The final step is to elute  $^{89}\text{Zr}$  with 1 M oxalic acid.

An elution distribution profile was determined with 10 x 200  $\mu\text{l}$  1 M oxalic acid injections into the module with a pump speed of  $\sim 3.5$  ml/min (Figure 3). To see any difference in separation trends for a smaller foil, the same procedure was also performed for irradiated 1.5 g foils. Elution with oxalic acid can be performed using either with a prefilled syringe connected to the module or via an external line out from the hot cell to the lab as in Figure 1. The latter gives the operator more flexibility if several [ $^{89}\text{Zr}$ ]Zr-oxalate batches are to be delivered.

## Labeling

The separated activity from the third 200  $\mu\text{l}$  fraction elute was used to label trastuzumab-Fab and tested against HER2 expressing SKOV-3 cells according to Siikanen et al.<sup>8</sup> The trastuzumab-Fab is obtained by fragmentation using a FragIT kit (Genovis, Lund).

## Specific activity

To calculate specific activity, the amount of cold zirconium was measured with inductively coupled plasma optical emission spectrometry (ICP-EOS) with a detection limit of 0.5 ppb for zirconium. For the ICP-EOS analysis the first 1 ml eluted oxalic acid separated from one un-irradiated foil was analysed.

## RESULTS

Separation of ten 200  $\mu\text{l}$  [ $^{89}\text{Zr}$ ]Zr-oxalate fractions took less than 90 min. More than 85% of the activity measured in foils was collected in total of 2 ml oxalic acid. No remarkable differences in elution profiles were observed for a 1.5g and a 3g foil separation. Approximately 90% of the activity transferred to the resin was eluted in the first milliliter of oxalic acid (Figure 3).

The overall labeling yield of  $^{89}\text{Zr}$ -Df-trastuzumab-Fab was over 99%. The radioconjugate was stable in PBS for at least 2 weeks and retained its binding specificity to HER2 receptors *in vitro* (Figure 4).

The ICP-EOS indicated a mass of 1.1  $\mu\text{g}$  of cold Zr.




FIGURE 3. Elution profile averaged over three separations




FIGURE 4. Binding specificity of  $^{89}\text{Zr}$ -trastuzumab-Fab on HER2-expressing SKOV-3 cells.  $^{89}\text{Zr}$ -trastuzumab-Fab binds specifically to HER2 receptors in SKOV-3 cells as shown by a pre-saturation of excess of cold trastuzumab (blocked) before addition of labeled substance

## DISCUSSION AND CONCLUSIONS

In this work a relatively large foil has been used, about 5-10 times more yttrium than common at other radiopharmaceutical production sites. Approximately 10 % of foil activity was trapped in the dissolving vial. This loss is expected to be lower for sites with smaller target masses which mean smaller dissolving volumes.

One hour irradiation yields about 2.1 GBq  $^{89}\text{Zr}$  in a 3g foil. Approximately 1.8 GBq (85%) of foil activity is eluted in 2 ml oxalic acid. About 1.6 GBq (90%) of total eluted activity is eluted in the first 1 ml oxalic acid fraction according to several separations like the example in Figure 3. Provided that irradiation conditions do not increase the amount of cold zirconium the specific activity of  $^{89}\text{Zr}$  is expected to be about 1.5 GBq/ $\mu\text{g}$  of zirconium (1 h irradiation). However more extensive studies involving specific activity and metal contamination measurements have been initiated to examine the separated product more thoroughly.

We have designed a simple and inexpensive automated separation module which also reduces personnel doses. This is helpful for sites which aim to distribute large amounts of  $^{89}\text{Zr}$  activity. The system is flexible and can be useful for other similar tasks if the setup and relay-pump sequences are reconfigured. Relay card, pump and pinch valves cost about 100 € each. All syringes, 3-way stopcocks and tubes in contact with liquids and foil are disposables.

## ACKNOWLEDGMENTS

The authors want to thank Jan Hultqvist from the workshop at Skåne University Hospital, Lund for machining the module parts. The authors also want to thank the John and Augusta Person Foundation, the Swedish Cancer Foundation, the Swedish Research Council, the Berta Kamprad Foundation, the Gunnar Nilsson Cancer Foundation and Medical faculty ALF grants for research support.

## REFERENCES

1. J. M. Link et al., *J. Labelled. Compd. Radiophar* **23**, 1297-1298 (1986).
2. G.A. van Dongen and M.J. Vosjan, *Cancer. Biother. Radiopharm* **25**, 375-385 (2010).
3. W.E. Meijs et al., *Appl Radiat Isot* **45**, 1143-1147 (1994).
4. G.W. Severin et al., *Med. Chem* **7**, 389-394 (2011).
5. M.J. Vojsan et al., *Nat. Protoc* **5**, 739-743 (2010).
6. J. Siikanen et al., "Multi-Purpose Solid Target System for a MC 17 Scanditronix Cyclotron" in *12<sup>th</sup> International Workshop on Targetry and Target Chemistry*, edited by J.M. Link, Seattle, WA, 2008, pp. 45-46.
7. J.P. Holland et al., *Nucl. Med. Biol* **36**, 729-739 (2009).
8. J. Siikanen et al., "Labeling proteins with  $^{89}\text{Zr}$  separated from large yttrium bulks" in *J Nucl Med Meeting Abstracts*, San Antonio, TX, 2011, **52** (Supplement 1):1504.