



Implementation of a Light Source in a TEM Sample Holder for In-situ Studies of Photocatalytic Materials

Cavalca, Filippo; Hansen, Thomas Willum; Wagner, Jakob Birkedal; Dunin-Borkowski, Rafal E.; Kardynal, Beata; Laursen, Anders Bo; Dionigi, Fabio; Dahl, Søren

Publication date:
2012

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):
Cavalca, F., Hansen, T. W., Wagner, J. B., Dunin-Borkowski, R. E., Kardynal, B., Laursen, A. B., Dionigi, F., & Dahl, S. (2012). *Implementation of a Light Source in a TEM Sample Holder for In-situ Studies of Photocatalytic Materials*. Abstract from 2011 MRS Spring Meeting and Exhibit, San Francisco, California, United States. <http://www.mrs.org/spring2011/>

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Implementation of a light source in a TEM sample holder for *in situ* studies of photocatalytic materials

F. Cavalca¹, T. W. Hansen¹, J. B. Wagner¹, R. E. Dunin-Borkowski¹, B. E. Kardynal², A. B. Laursen³, F. Dionigi³, S. Dahl³

¹ Center for Electron Nanoscopy, Technical University of Denmark, DK-2800 Lyngby, Denmark

² Department of Photonics Engineering, Technical University of Denmark, Ørstedes Plads 343, DK-2800 Lyngby, Denmark

³ Center for Individual Nanoparticle Functionality, Technical University of Denmark, Fysikvej, Building 312, DK-2800 Lyngby, Denmark

Photocatalysts are of fundamental interest for sustainable energy research [1]. By means of transmission electron microscopy (TEM) it is possible to obtain insight into the structure, composition and reactivity of photocatalysts. Such insight can be used for their further optimization [2].

We have constructed a specimen holder capable of shining light onto samples inside the TEM. The holder contains a laser diode and an optical system that guides light onto a sample with maximum power transmission. The source can be changed and tuned, in principle spanning the whole visible and UV spectrum. The device can be used inside an environmental TEM (ETEM) allowing specimens to be analyzed during exposure to a controlled gas atmosphere and illumination.

The holder is presently being used to study a variety of photoreactive materials and structures, including photocatalysts, photonic devices and solar cells. For example, electron holography is being used to study p-n junctions both in the presence and in the absence of light in order to assess electron beam induced charging and discharging effects during laser light exposure [3].

Here, we present results from ETEM studies of light-induced phenomena that include metal nanoparticle photodeposition, light-driven particle discharging and photodegradation. We concentrate on phase transitions of Cu₂O nanocubes under visible light exposure in the presence of water vapor, which we study *in situ*. Cu₂O is an active photocatalyst for water splitting under visible light irradiation, but it undergoes photodegradation in an aqueous environment [4].

[1] Herrmann, J. M., *Top. Catal.* **2005**, 34, (1-4), 49-65.

[2] Tsujimoto, M., S. Moriguchi, et al., *J. Electron. Microsc. (Tokyo)* **1999**, 48, (4), 361-366.

[3] Shindo, D., K. Takahashi, et al., *J Electron Microsc. (Tokyo)* **2009**, 58, (4), 245-249.

[4] Nagasubramanian, G., A. S. Gioda, *J. Electrochem. Soc.* **1981**, 128, (10), 2158-2164.