



Speciation analysis of trace elements in foodstuffs by HPLC- and GC-ICPMS

Sloth, Jens Jørgen

Publication date:
2012

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):
Sloth, J. J. (Author). (2012). Speciation analysis of trace elements in foodstuffs by HPLC- and GC-ICPMS. Sound/Visual production (digital)

General rights


Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.


- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Speciation analysis of trace elements in foodstuffs by HPLC- and GC-ICPMS

Jens J. Sloth





4 metals on the WHO IPCS list of top-10 priority chemicals of public concern

CURRENT SITUATION IN EU LEGISLATION:


Foodstuffs

MLs for Pb, Cd, Hg and Sn
EU directive 2006/1881/EC

Animal feedingstuffs

MLs for As, Pb, Cd and Hg
EU directive 2002/32/EC

Only maximum levels for
TOTAL CONCENTRATION
of the metals



Analysis of total content vs. speciation analysis

~~Analysis of total concentration~~

Arsenic
Tin
Mercury
Lead
Chromium
Selenium
Iodine


Aim: get information on e.g.:

- biological activity
- toxicity
- mc
- bic
- lifetime, fate and metabolism
- chemical and physical activity

these characteristics are species related


Speciation analysis is required!!

Speciation analysis – chromatography and element-specific detection




EXAMPLES

- Arsenic speciation analysis by HPLC-ICPMS
- Organotin speciation analysis by GC-ICPMS
- Selenium speciation by HPLC-ICPMS
- Nanoparticle determination by FFF-MALS-ICPMS





- Beyond total element determination
- Species-specific information




Arsenic compounds in the marine environment


More than **50** different arsenic species have been found in the marine environment




Dimethylarsinoylriboside




Dimethyl As-sugar 1






Dimethyl As-sugar 2



Dimethyl As-sugar 3



Dimethyl As-sugar 4



Trimethylarsonioriboside


R = 1,2,3,4,5

Speciation analysis of arsenic of scallop kidney

Cation-exchange with gradient elution – extraction with aqueous methanol


Column: Chrompack Ionospher 5C; Mobile phase: Pyridine; pH = 2.7

- seven compounds identified by coelution with available standards
- **16** non-identified peaks



Example – arsenic speciation

Important for correct risk assessment



~ 40 – 250 μg inorg As

1 kg rice => 50-300 μg As



<10 μg inorg As

1 kg fish => 3000-10000 μg As


There is most focus on rice from a food safety point of view – why???



The chemical form of arsenic is important
=> Arsenic speciation


Focus
on
Food
Safety

Arsenic compounds in the marine environment




**Most
toxic
form of
arsenic!**


Inorganic arsenic



Arsenous acid
As(III)



Arsenic acid
As(V)



EFSA (2009) and JECFA (2010) opinions on arsenic in food

- Old PTWI value (WHO, 1988) was withdrawn ($15 \mu\text{g/kg bw/week}$)




- **NEW!** $\text{BMDL}_{1.0} = 0.3 - 8 \mu\text{g/kg bw per day}$ for inorganic arsenic
- \Rightarrow EU dietary exposures within this range for average and high level consumers
- \Rightarrow Risk to some consumers cannot be excluded



- **NEW!** $\text{BMDL}_{0.5} = 3 \mu\text{g/kg bw per day}$ for inorganic arsenic
 $\Rightarrow 0.5\%$ increased incidence of lung cancer for 12 y exposure

- “...there is a need to produce speciation data for different food commodities to support dietary exposure assessment...”
- “...more accurate information on the inorganic arsenic content of foods is needed to improve assessments of dietary exposures to inorganic arsenic”
- “...need for validated methods for selective determination of inorganic arsenic in food matrices”

Inorganic arsenic in wild caught fish => no concern



Norwegian survey

900 individual fish samples

- Atlantic halibut
- Cod
- Greenland halibut
- Mackerel
- Herring
- Tusk

Results


Total arsenic.....0.3-110 mg/kg

Inorganic arsenic.... < 0.01 mg/kg

(only 37 samples > LOQ)



...but in bivalves high contents in some samples...



Data from 175 blue mussel (*Mytilus edulis*) samples collected along the Norwegian Coastline.



Arsenic in rice – an emerging health issue?

Environmental Pollution 152 (2008) 746–749

Rapid communication

Inorganic arsenic levels in baby rice are of concern

Andrew A. Meharg^{a,*}, Guoxin Sun^b, Paul N. Williams^{a,b}, Eureka Adomako^a,
Claire Deacon^a, Yong-Guan Zhu^b, Joerg Feldmann^c, Andrea Raab^c

^a School of Biological Sciences, University of Aberdeen, Cruickshank Building, St. Machar Drive, Aberdeen AB24 3UU, UK


^b Research Centre for Eco-Environmental Sciences, Chinese Academy of Sciences, Beijing 100085, China

Median consumption of organic arsenic levels for UK babies from baby rice is above threshold considered safe.

- 17 samples from supermarkets in Aberdeen
- Total arsenic levels: 0.12 – 0.47 mg/kg
- Inorganic arsenic: 0.06 – 0.16 mg/kg (33 – 69 % of tAs)
- 35% above Chinese max level of 0.15 mg/kg iAs
- No regulation on As in food in EU (yet!)



Arsenic in rice products DK - survey




105 samples

- white rice
- brown rice
- red rice
- black rice
- rice crackers

33 samples > 0,15 mg/kg

- 2 parboiled (20%)
- 4 brown (50%)
- 4 red (50%)
- 5 black (71%)
- 1 Basmati (10%)
- 1 Pudding rice (9%)
- 1 wild rice (20%)
- 15 rice crackers (100%)

Arsenic in rice products DK - survey



105 samples

- white rice
- brown rice
- red rice
- black rice
- rice crackers

22 samples > 0,2 mg/kg

- 1 parboiled (10%)
- 1 brown (12%)
- 3 red (37%)
- 2 black (28%)
- 0 Basmati (0%)
- 0 Pudding rice (0%)
- 0 wild rice (0%)
- 15 rice crackers

Rice cracker mean: 0.31 mg/kg – intake 50 g/dag => 15 µg iAs (~1 µg/kg bw/dag)

> EFSA BMDL₀₁ 0.3-8 µg/kg bw/dag



Inorganic arsenic in chinese food supplements

Name of Food supplement	Total Arsenic ($\mu\text{g/g}$)	Inorganic arsenic ($\mu\text{g/g}$)
Xiao Yao Wan	0.82	0.85
Bu Zhong Yi Qi Wan	0.62	0.50
Da Bu Yin Wan	0.59	0.55
Six Flavor teapills	0.72	N.D.
Golden Book Teapills	0.58	0.57
Xiang Sha Liu Jun Zi Wan	0.94	0.80
Gan Mao Ling	1.24	1.01
Chuan Xin Lian	5.00	3.17
Bi Yan Pian	0.70	0.58
Arouse power	1.12	1.02
Bio Chlorella	0.62	0.21
Unik Spirulina Kapsler	2.59	0.13
Chlorella	0.58	0.03
Ez-Biloba	0.63	0.67
Qvinde Dong Quai	0.68	0.48




Gan Mao Ling
(against flu and common cold)

Rec dose: 18 pills per day
 \Rightarrow iAs $\sim 13 \mu\text{g/day}$
 $\Rightarrow 0.22 \mu\text{g/kg bw/day}$ (@60 kg)
 Close to EFSA BMDL_{01} !!


Arsenolipids – a class of lipidsoluble compounds

Sand eel oil (important in fishfeed production)

GC-ICPMS




Arsenic hydrocarbons



- Still limited knowledge about:
 - chemical structures
 - biogenesis
 - bioavailability, metabolism
 - toxicity??




Uses of organotin compounds 1 - pesticides



- **Triphenyltin acetate or hydroxide** (Brestan, Aquatin or Telostan) e.g. used for control the population of brackish water snails *Cerithidea cingulata* in the pond culture of milkfish



- **Acaricide azocyclotin**
- used against spiders, mites on e.g. be
- **Fenbutatin oxide**
- used as insecticide on almonds, grape



Uses of organotin compounds 2 - antifoulants



- Prevention of growth of plants and animals on ships and stationary marine constructions
- 1991: Banned on yachts <25 m in the Nordic countries
- 1999: Banned on yachts in all EU
- 2003: TBT ban on re-painting bigger ships with TBT paint
- 2008: TBT paint must be sealed or not used at all in EU ships




Example of metal plate treated with both antifouling paint and non-antifouling paint (www.ortepa.org)

Routes of exposure to organotin compounds

Used in

- Agriculture
- Antifoulings
- Industry



Data from seafood surveys at DTU Food

- organotin compounds in marine samples

- **Project on various fishtypes:**
 - 214 samples of fish and shellfish from Danish fish wholesales
 - Generally sub-ppb levels found for all OT
 - Herring (up to 19 μg OT/kg)
 - Shark (up to 39 μg OT/kg)
 - Eel (up to 3200 μg OT/kg) !!
- **Projects on bivalves from Denmark:**
 - mono- , di- and tri-butylated and phenylated compounds
 - 100 bivalves samples/yr
 - generally only low levels found (sum OT < 10 μg /kg)



Routes of exposure to organotin compounds

Used in


- Agriculture
- Antifoulings
- Industry

- PVC-Materials

TDI: 0.25 $\mu\text{g}/\text{kg bw}/\text{day}$

$\sum \text{TBT, DBT, TPhT and DOT}$

$\Rightarrow \sim 17 \mu\text{g}/\text{day} @ \text{bw } 70\text{kg}$



Uses of organotin compounds 3 – PVC stabilisers

- PVC is a versatile polymer used worldwide for a variety of products
- At processing temperature PVC is unstable
- Stabilisers are added – functions against heat and light
- Mono- or di-alkyl tin compounds
- RSnX_3 and R_2SnX_2
- R = methyl-, butyl-, octyl-, dodecyl-
- X = mercaptoester, carboxylate, sulfide

pipes



containers



windows



cling films



Food Contact Materials



Applications of rigid and flexible PVC with OT

Rigid Applications	Quantity (t/yr)
Packaging, incl food contact materials , credit cards	12343
Rigid construction incl foamed sheeting	1016
Thin rigid film	290
Bottles	290
Pipes and moulding	290
Profile extrusions (e.g. windows)	290
Flexible Applications	
Flooring	312
Wallcovering	312
Steel coating	312
Miscellaneous (e.g. T-shirt printing)	156

Legislation on OTCs in Food Contact Materials




Compounds	Maximum level ($\mu\text{g Sn/kg foodstuff}$)
\sum DBT, TBT, TPhT and DOT	40 (6)
\sum MMT, DMT	180
MOT	1200
MDDT	12000 (50)
DDDT	24000 (50)

Ref: EFSA (2005); proposed EFSA values in parenthesis

Assumptions:

- 1 kg food per 6 dm²
- 100 mL in contact with 0.6 dm²

- Max levels on organotin migrating from the packaging material
- Testing by the use of food simulators (water, acid, oil, alcohol etc)
- BUT no maximum levels on organotins in the foodstuff itself!!




Organotin migration from Food Contact Materials II

Small scale survey on 33 FCMs

Baking paper, PVC cling films, silicone baking forms, lids with PVC gaskets

PUR-agglomerated cork wine stoppers



- PVC lid
- 3% acetic acid

Overlaid standard and sample

➤ DBT concentration: 9.9 µg/kg

TIC: 005SMPL.D
TIC: 014SMPL.D

DBT

> EFSA guideline
value of 6 µg/kg

Output of DK survey:

- 33 samples
- 11 contained OT (mainly DBT)
- 3 exceeded EFSA guideline limit


Selenium in commercial food supplements

Organic bound Se?



Selenite?

Selenate?



Se yeast?


Selenomethionine?

Amino chelated Se?

Is the selenium source declared correctly??

Se speciation by HPLC-ICPMS




Supplement declared as organic bound Selenium (125 µg/tablet)



Silver nanoparticles in products related to food and beverages



The analytical platform



**Asymmetric flow
field flow
fractionation**

Optical detection
(multi angle and
dynamic light
scattering, UV and
fluorescence)

**Inductively
coupled plasma
mass spectrometry
(ICP-MS)**

Particle separation
according to their size
(small NPs elute first)

Particle detection
(fractogram)



Size determination
(root mean square,
hydrodynamic and
geometric radius)

Elemental detection
for identification of
particles

Quantification

Determination of nanoparticle size distribution

Fraction collection for transmission electron microscopy (TEM)



Perspective – trace elements in nanoform

NanoLyse Project "Nanoparticles in Food: Analytical methods for detection and characterisation"

Validated methods for the determination of inorganic ENP in food extracts, based on size separation, size determination and specific detection

silver nanoparticles in lean meat



silica nanoparticles in tomato soup



Acknowledgements and funding sources

Coworkers:

Rie R. Rasmussen, Rikke V. Hedegaard, Bjørn Schmidt,
Xenia T. Trier, Katrin Löschner and Erik H. Larsen

Inge Rokkjær, Gudrun Hilbert and Dorte L. Cederberg

Kåre Julshamn and A.K. Lundebye



Funding sources:


European Community's Seventh Framework Programme

Danish Food Administration (DFVA)

Norwegian Research Council



Thanks for your attention! 



**6th Nordic Conference
on
Plasma Spectrochemistry
June 10-13, 2012, Loen, Norway**

www.nordicplasma.com

