



Early Detection Of Chronic Bacterial Infections By Microsensors

Al Atraktchi, Fatima Al-Zahraa; Molin, Søren; Johansen, Helle Krogh; Svendsen, Winnie Edith

Publication date:
2015

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):

Al Atraktchi, F. A-Z., Molin, S., Johansen, H. K., & Svendsen, W. E. (2015). *Early Detection Of Chronic Bacterial Infections By Microsensors*. Poster session presented at The 3rd International Conference on Nanotechnology in Medicine (NANOMED), Manchester, United Kingdom.

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

EARLY DETECTION OF CHRONIC BACTERIAL INFECTIONS BY MICROSENSORS



Fatima AlZahra'a Alatraktchi^{1,2,3}, Søren Molin^{2,3}, Helle Krogh Johansen^{3,4}, Winnie E. Svendsen¹

¹ Department of Micro- and Nanotechnology - Technical University of Denmark

² Department of Systems Biology - Technical University of Denmark


³ Novo Nordisk Foundation Center for Biosustainability - Technical University of Denmark

⁴ Cystisk Fibrose Klinikken & Klinisk Mikrobiologisk Afdeling - Rigshospitalet

1 Motivation

The bacterial pathogen *Pseudomonas aeruginosa* is involved in many dangerous infections. A precise and early diagnosis of an infection is essential for successful eradication treatment. Thus, there is need for fast and sensitive methods, which can be used in the clinic to diagnose bacterial infections in patients.


The redox-active pyocyanin is produced by *P. aeruginosa* prior to virulence activity. Monitoring of pyocyanin as an infection biomarker could consequently enable early detection of bacterial infections. This study has shown that electrochemical sensors can detect pyocyanin levels in the nanomolar range, which will enable early diagnosis of the bacteria before chronic establishment.




P. aeruginosa produce pyocyanin prior to virulent activity. Ref: Free stock

2 Concept


Micro-Electro-Mechanical Systems (MEMS) is a technology with devices ranging from micrometer to several millimeter scale. Commercial electrochemical microsensors were used to detect pyocyanin by cyclic voltammetry, chronoamperometry and square wave voltammetry using a three-electrode configuration.



Cyclic voltammetry of pyocyanin




Chronoamperometry of pyocyanin




Square Wave Voltammogram of pyocyanin

3 Status

It is a goal to directly detect pyocyanin in sputum samples from patients to diagnose the stage of infection. Identification of pyocyanin in clinically relevant media like lysogeny broth (LB) used for bacterial culturing is essential. Detection and quantification of pyocyanin has been conducted in both LB and Artificial Sputum (ASM) using square wave voltammetry.



Square Wave Voltammogram of pyocyanin in LB medium (blue) compared to pure LB medium (black)



Square Wave Voltammogram of pyocyanin in ASM (blue) compared to pure ASM (black)

4 Impact

This study is expected to reveal valuable information regarding the progress of an infection. The technique will assist medical doctors in designing individual antibiotic treatment plans for each specific situation. An enhanced antibiotic treatment strategy could contribute to a longer lifetime of patients.

Pyocyanin is expected to be detected as a label-free biomarker for *P. aeruginosa* in point-of-care systems in the clinic or for home use to monitor the condition of the patients.



Mobile application to display measurements of pyocyanin. Ref: Francois Patou, DTU