



## Sculpted light for twophoton optogenetics using GPC

Glückstad, Jesper

*Publication date:*  
2013

*Document Version*  
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

*Citation (APA):*  
Glückstad, J. (2013). *Sculpted light for twophoton optogenetics using GPC.*

---

### General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

# Sculpted light for two-photon optogenetics using GPC

J. Glückstad  
DTU Fotonik, Dept. of Photonics Engineering,  
Programmable Phase Optics  
Technical University of Denmark,  
Ørsted Plads 343, 032  
DK-2800 Kgs. Lyngby, Denmark

[jesper.gluckstad@fotonik.dtu.dk](mailto:jesper.gluckstad@fotonik.dtu.dk)

[www.ppo.dk](http://www.ppo.dk)

The synergy between nanotechnology, biotechnology and photonics is spawning the emerging fields of nano-biotechnology and nano-biophotonics. Photonic innovations already hurdle the diffraction barrier for imaging with nanoscopic resolutions. However, scientific hypothesis testing demands tools, not only for observing nanoscopic phenomena, but also for reaching into and manipulating nanoscale constituents in this domain. This talk is describing the new use of proprietary strongholds we are currently establishing on new means of sculpting light for bio-probing and excitation at the smallest cellular scales.

Two-photon spatio-temporal light engineering can be used to expand the microscopic imaging modalities available to assist this light-driven nano-manipulation approach. Featured in a recent issue of Nature Methods we were pioneering research in neuro-photonics and optogenetics highly useful for future biophotonics undertakings on the smallest cellular scales. This research promises a powerful approach for controlling light-gated ion channels and pumps that facilitates the probing of intact neural circuits by manipulating the activity of groups of genetically similar neurons. This makes it possible to precisely aim space-time sculpted light at single neuronal processes, neurons or groups of neurons. The underlying light-engineering is currently being scientifically refined in a strong international context to address - for the first time - arbitrary and speckle-free '4D' spatial and temporal two-photon light-sculpting.

## References

1. P. Rodrigo, L. Kelemen, D. Palima, C. Alonzo, P. Ormos, J. Glückstad, "Optical microassembly platform for constructing reconfigurable microenvironments for biomedical studies," *Opt. Express* 17, 6578-6583 (2009).
2. J. Glückstad, "Sorting particles with light," *Nature Materials* 3, 9-10 (2004).
3. J. Glückstad, "Optical manipulation: Sculpting the object," *Nature Photonics* 5, 7-8 (2011).
4. E. Papagiakoumou, F. Anselmi, A. Begue, V. de Sars, J. Glückstad, E. Isacoff, V. Emiliani, "Scanless two-photon excitation of channelrhodopsin-2," *Nature Methods* 7, 848-854 (2010).
5. D. Palima, A. R. Bañas, G. Vizsnyiczai, L. Kelemen, P. Ormos, and J. Glückstad, "Wave-guided optical waveguides," *Opt. Express* 20, 2004-2014 (2012).
6. D. Palima and J. Glückstad, "Gearing up for optical microrobotic manipulation: mechanical actuation of synthetic microstructures by optical forces," *Laser & Photonics Reviews* 7, 478-494 (2013).
7. H. Ulriksen, J. Thøgersen, S. Keiding, I. Perch-Nielsen, J. Dam, D. Palima, H. Stapelfeldt, and J. Glückstad, "Independent trapping, manipulation and characterization by an all-optical biophotonics workstation," *J. Eur. Opt. Soc-Rapid* 3, 08034 (2008).